


Auto-Smear™ Microscope Slides


Recommended for automatic blood smear, including Sysmex fully automated slide-maker Sp1000i, COULTER® SlideMaker & SlideStainer Hematology System LH 755, UniCel® DxH SlideMaker & SlideStainer Coulter Cellular Analysis System and COULTER® Hematology Automation

*Markable with the conventional labeling system, pencil, mark pens as well as the popular 2-D barcode based slide printers, such as the Inkjet IPS (Leica), Slide-Mate® (Thermo Fisher) and Signature® (Primer) slide printers



Dimensions: 25.0mmx75.0mm
Thickness: 1.1±0.1mm

Cleanliness and wettability of slide surfaces are strictly and uniformly controlled



- Color frosted area of approx. 20 mm depth, displays at one surface end
- Color frosted end is resistant to all common chemicals and routine stains that are utilized in the laboratory

CSW

MARKING AREA	P/N	PACKING
PINK	5860-0001	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
PINK	5860-0002	72pcs/box, 1,440pcs/carton
WHITE	5860-0011	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
WHITE	5860-0012	72pcs/box, 1,440pcs/carton
GREEN	5860-0045	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
GREEN	5860-0046	72pcs/box, 1,440pcs/carton

SUPER GRADE

MARKING AREA	P/N	PACKING
PINK	5860-0101	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
PINK	5860-0102	72pcs/box, 1,440pcs/carton
WHITE	5860-0111	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
WHITE	5860-0112	72pcs/box, 1,440pcs/carton
GREEN	5860-0145	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
GREEN	5860-0146	72pcs/box, 1,440pcs/carton

PapSmear™ Microscope Slides

Suitable for manual blood smear preparations, also ideal raw material for liquid-based cytology, e.g., TCT & LCT slide preparation

- Dimensions: 25.0mmx75.0mm, Thickness: 1.1±0.1mm
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requirements

CSW

-50/PK 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	Ground Edges 45°corners	Beveled Edges 45°corners
WHITE	5860-0015	5860-0017

-72/PK 72pcs/box, 1,440pcs/carton

MARKING AREA	Ground Edges 45°corners	Beveled Edges 45°corners
WHITE	5860-0016	5860-0018

SUPER GRADE

-50/PK 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	Ground Edges 45°corners	Beveled Edges 45°corners
WHITE	5860-0115	5860-0117

-72/PK 72pcs/box, 1,440pcs/carton

MARKING AREA	Ground Edges 45°corners	Beveled Edges 45°corners
WHITE	5860-0116	5860-0118

Instruction for Blood Smear

A Place a drop on the end of slide, Hold the spatula at a 30° ~45°angle, toward the drop of blood

B Wait until the blood spreads along the entire width of the spatula

C Hold the spatula at the same angle, push it forward rapidly and smoothly

