

Adhesion Microscope Slides

The continuous liquid-bath in stains as well as chemicals and the frequent washing throughout the immunohistochemical staining workflow, especially the immediate washing after the antigen-retrieval heating do cause great risk of tissue section loss or cell loss. The surface of the **AA**® Color-PLUS™ microscope slides is coated with a permanent especially-formulated adhesion which reliably binds the tissue sections to the glass, considerably minimizing the possibility of tissue-section or cell loss.

The **AA**® Color-PLUS™ microscope slides are also highly recommended in routine H&E staining to secure the problematic tissue sections like brain and spinal cord specimens.

CAUTIONS

- 1.Store the products under dry and room temperature conditions.
- 2. Use the products as soon as possible when the package is open.
- 3.The color frosted layer may fall off the glass where caustic conditions are applied.
- 4. Conduct a trail prior to use the products for the first time.

Applications

Recommended ★
Strongly recommended ★ ★

PCB

PCL

PCC

PCI

HISTOLOGY AND PATHOLOGY

Paraffin embedded section
Frozen section
Liquid-rich specimens



★ ★
★ ★
★

★ ★
★ ★
★

★ ★
★ ★
★

★ ★
★ ★
★ ★

AUTOMATED IHC STAINER

Roche Ventana BenchMark
Leica BOND-MAX
Dako Autostainer Link 48

★ ★
★ ★
★ ★

★ ★
★ ★
★ ★

★ ★
★ ★
★ ★

★ ★
★ ★
★ ★

HIER PROCESSES

Acid Retrieval Solution
Alkaline Retrieval Solution

★ ★
★ ★

★ ★
★ ★

★ ★
★ ★

★ ★
★ ★

CYTOLOGY

LBC, membrane
LBC, centrifugation
Microbiology
Hematology

★ ★
★ ★
★ ★
★

★ ★
★ ★
★ ★
★

★ ★
★ ★
★ ★
★

★ ★
★ ★
★ ★
★




Color-Plus™ PCB Microscope Slides

Recommended for manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil and mark pens, but not supportive for direct printing with the 2-D barcode based printers.

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophobic surface promotes the adhesion of specimens
- Color frosted area of approx. 20 mm depth, displays at one surface end
- Color frosted end is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requirements



-50/PK
Packaging: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	 Ground Edges 90°corners		 Ground Edges 45°corners		 Beveled Edges 45°corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7141-02	0303-7141-02	0312-7141-02	0302-7141-02	0314-7141-02	0304-7141-02
Green	0313-7141-08	0303-7141-08	0312-7141-08	0302-7141-08	0314-7141-08	0304-7141-08
Orange	0313-7141-10	0303-7141-10	0312-7141-10	0302-7141-10	0314-7141-10	0304-7141-10
Pink	0313-7141-12	0303-7141-12	0312-7141-12	0302-7141-12	0314-7141-12	0304-7141-12
White	0313-7141-16	0303-7141-16	0312-7141-16	0302-7141-16	0314-7141-16	0304-7141-16
Yellow	0313-7141-17	0303-7141-17	0312-7141-17	0302-7141-17	0314-7141-17	0304-7141-17

-72/PK
Packaging: 72pcs/box, 1,440pcs/carton

MARKING AREA	 Ground Edges 90°corners		 Ground Edges 45°corners		 Beveled Edges 45°corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7151-02	0303-7151-02	0312-7151-02	0302-7151-02	0314-7151-02	0304-7151-02
Green	0313-7151-08	0303-7151-08	0312-7151-08	0302-7151-08	0314-7151-08	0304-7151-08
Orange	0313-7151-10	0303-7151-10	0312-7151-10	0302-7151-10	0314-7151-10	0304-7151-10
Pink	0313-7151-12	0303-7151-12	0312-7151-12	0302-7151-12	0314-7151-12	0304-7151-12
White	0313-7151-16	0303-7151-16	0312-7151-16	0302-7151-16	0314-7151-16	0304-7151-16
Yellow	0313-7151-17	0303-7151-17	0312-7151-17	0302-7151-17	0314-7151-17	0304-7151-17

PCL Color-Plus™ PCL Microscope Slides **NEW**

Recommended for manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil, mark pens as well as the 2-D barcode based printers

Adhesion Microscope Slides for Printing Printers

Recommended for marking with the well-known brand printers, including the Inkjet IPS (Leica), Slide-Mate® (Thermo Fisher) and Signature® (Primer) slide printers

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophobic surface promotes the adhesion of specimens
- Color frosted area of approx. 20 mm depth, displays at one surface end
- Color frosted end is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requirements



-50/PK
Packaging: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
 Blue	0313-7181-02	0303-7181-02	0312-7181-02	0302-7181-02	0314-7181-02	0304-7181-02
 Green	0313-7181-08	0303-7181-08	0312-7181-08	0302-7181-08	0314-7181-08	0304-7181-08
 Orange	0313-7181-10	0303-7181-10	0312-7181-10	0302-7181-10	0314-7181-10	0304-7181-10
 Pink	0313-7181-12	0303-7181-12	0312-7181-12	0302-7181-12	0314-7181-12	0304-7181-12
 White	0313-7181-16	0303-7181-16	0312-7181-16	0302-7181-16	0314-7181-16	0304-7181-16
 Yellow	0313-7181-17	0303-7181-17	0312-7181-17	0302-7181-17	0314-7181-17	0304-7181-17

-72/PK
Packaging: 72pcs/box, 1,440pcs/carton

MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
 Blue	0313-7191-02	0303-7191-02	0312-7191-02	0302-7191-02	0314-7191-02	0304-7191-02
 Green	0313-7191-08	0303-7191-08	0312-7191-08	0302-7191-08	0314-7191-08	0304-7191-08
 Orange	0313-7191-10	0303-7191-10	0312-7191-10	0302-7191-10	0314-7191-10	0304-7191-10
 Pink	0313-7191-12	0303-7191-12	0312-7191-12	0302-7191-12	0314-7191-12	0304-7191-12
 White	0313-7191-16	0303-7191-16	0312-7191-16	0302-7191-16	0314-7191-16	0304-7191-16
 Yellow	0313-7191-17	0303-7191-17	0312-7191-17	0302-7191-17	0314-7191-17	0304-7191-17

PCL Color-Plus™ PCL Microscope Slides **NEW**


Recommended for manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil, mark pens as well as the 2-D barcode based printers

Adhesion Microscope Slides for Etching Printers


Recommended for marking with the slide etching labelers

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophobic surface promotes the adhesion of specimens
- Color frosted areas of approx. 20 mm depth, distribute on both sides but at the same slide end for marking
- A black color opposite the color side makes the marking clearer
- Color frosted areas are resistant to all common chemicals and routine stains that are utilized in the laboratory
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requests and experiments




-72/PK
Packaging: 72pcs/plastic box, 1,440pcs/carton


MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
 Blue	0313-3291-02	0303-3291-02	0312-3291-02	0302-3291-02	0314-3291-02	0304-3291-02
 Green	0313-3291-08	0303-3291-08	0312-3291-08	0302-3291-08	0314-3291-08	0304-3291-08
 Orange	0313-3291-10	0303-3291-10	0312-3291-10	0302-3291-10	0314-3291-10	0304-3291-10
 Pink	0313-3291-12	0303-3291-12	0312-3291-12	0302-3291-12	0314-3291-12	0304-3291-12
 White	0313-3291-16	0303-3291-16	0312-3291-16	0302-3291-16	0314-3291-16	0304-3291-16
 Yellow	0313-3291-17	0303-3291-17	0312-3291-17	0302-3291-17	0314-3291-17	0304-3291-17




Paraffin Embedded Sections




Frozen Sections




IHC



Human Blood Cells



Cervical Samples



Yeast Cells

PCI Color-Plus™ PCI Microscope Slides



Recommended for automatic IHC staining with the Roche Ventana IHC automated Stainer (priviledged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil, mark pens as well as the 2-D barcode based slide printers

Adhension Microscope Slides for Printing Printers

Recommended for marking with the well-known brand printers, including the Inkjet IPS (Leica), Slide-Mate® (Thermo Fisher) and Signature® (Primera) slide printers

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophilic surface promotes the adhesion of specimens
- Color frosted area of approx. 20 mm depth, displays at one surface end
- Color frosted end is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requests and experiments



-50/PK

Packaging: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7161-02	0303-7161-02	0312-7161-02	0302-7161-02	0314-7161-02	0304-7161-02
Green	0313-7161-08	0303-7161-08	0312-7161-08	0302-7161-08	0314-7161-08	0304-7161-08
Orange	0313-7161-10	0303-7161-10	0312-7161-10	0302-7161-10	0314-7161-10	0304-7161-10
Pink	0313-7161-12	0303-7161-12	0312-7161-12	0302-7161-12	0314-7161-12	0304-7161-12
White	0313-7161-16	0303-7161-16	0312-7161-16	0302-7161-16	0314-7161-16	0304-7161-16
Yellow	0313-7161-17	0303-7161-17	0312-7161-17	0302-7161-17	0314-7161-17	0304-7161-17

-72/PK

Packaging: 72pcs/box, 1,440pcs/carton

MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7171-02	0303-7171-02	0312-7171-02	0302-7171-02	0314-7171-02	0304-7171-02
Green	0313-7171-08	0303-7171-08	0312-7171-08	0302-7171-08	0314-7171-08	0304-7171-08
Orange	0313-7171-10	0303-7171-10	0312-7171-10	0302-7171-10	0314-7171-10	0304-7171-10
Pink	0313-7171-12	0303-7171-12	0312-7171-12	0302-7171-12	0314-7171-12	0304-7171-12
White	0313-7171-16	0303-7171-16	0312-7171-16	0302-7171-16	0314-7171-16	0304-7171-16
Yellow	0313-7171-17	0303-7171-17	0312-7171-17	0302-7171-17	0314-7171-17	0304-7171-17

PCI Color-Plus™ PCI Microscope Slides




Recommended for automatic IHC staining with the Roche Ventana IHC automated Stainer (priviledged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil, mark pens as well as the 2-D barcode based slide printers

Adhension Microscope Slides for Etching Printers

Recommended for marking with the slide etching labelers


- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophilic surface promotes the adhesion of specimens
- Color frosted areas of approx. 20 mm depth, distribute on both sides but at the same slide end for marking
- Color frosted areas are resistant to all common chemicals and routine stains that are utilized in the laboratory
- A black color opposite the color side makes the marking clearer
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requests and experiments




-72/PK

Packaging: 72pcs/plastic box, 1,440pcs/carton


MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3271-02	0303-3271-02	0312-3271-02	0302-3271-02	0314-3271-02	0304-3271-02
Green	0313-3271-08	0303-3271-08	0312-3271-08	0302-3271-08	0314-3271-08	0304-3271-08
Orange	0313-3271-10	0303-3271-10	0312-3271-10	0302-3271-10	0314-3271-10	0304-3271-10
Pink	0313-3271-12	0303-3271-12	0312-3271-12	0302-3271-12	0314-3271-12	0304-3271-12
White	0313-3271-16	0303-3271-16	0312-3271-16	0302-3271-16	0314-3271-16	0304-3271-16
Yellow	0313-3271-17	0303-3271-17	0312-3271-17	0302-3271-17	0314-3271-17	0304-3271-17




Paraffin Embedded Sections




Frozen Sections




IHC



Human Blood Cells



Cervical Samples



Yeast Cells

PCC Color-Plus™ PCC Microscope Slides

Recommended for routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil, mark pens as well as the 2-D barcode based slide printers

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophilic surface promotes the adhesion of specimens
- Color frosted area of approx. 20 mm depth, displays at one surface end
- Color frosted end is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requests and experiments



-50/PK
Packaging: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
 Blue	0313-7183-02	0303-7183-02	0312-7183-02	0302-7183-02	0314-7183-02	0304-7183-02
 Green	0313-7183-08	0303-7183-08	0312-7183-08	0302-7183-08	0314-7183-08	0304-7183-08
 Orange	0313-7183-10	0303-7183-10	0312-7183-10	0302-7183-10	0314-7183-10	0304-7183-10
 Pink	0313-7183-12	0303-7183-12	0312-7183-12	0302-7183-12	0314-7183-12	0304-7183-12
 White	0313-7183-16	0303-7183-16	0312-7183-16	0302-7183-16	0314-7183-16	0304-7183-16
 Yellow	0313-7183-17	0303-7183-17	0312-7183-17	0302-7183-17	0314-7183-17	0304-7183-17


-72/PK
Packaging: 72pcs/box, 1,440pcs/carton

MARKING AREA	 Ground Edges 90° corners		 Ground Edges 45° corners		 Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
 Blue	0313-7193-02	0303-7193-02	0312-7193-02	0302-7193-02	0314-7193-02	0304-7193-02
 Green	0313-7193-08	0303-7193-08	0312-7193-08	0302-7193-08	0314-7193-08	0304-7193-08
 Orange	0313-7193-10	0303-7193-10	0312-7193-10	0302-7193-10	0314-7193-10	0304-7193-10
 Pink	0313-7193-12	0303-7193-12	0312-7193-12	0302-7193-12	0314-7193-12	0304-7193-12
 White	0313-7193-16	0303-7193-16	0312-7193-16	0302-7193-16	0314-7193-16	0304-7193-16
 Yellow	0313-7193-17	0303-7193-17	0312-7193-17	0302-7193-17	0314-7193-17	0304-7193-17

Liquid-Based Cytology Microscope Slides

Supplied for liquid-based cytology, e.g., TCT & LCT slide preparation.

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Coated and non-coated available
- Ground Edges, 45°clipped corners
- Frosted area of approx. 20mm depth, displays at one surface end
- Printed Circle or Bracket optional



Microscope Slides with Printed Circles

Compatible with single or double Cellfunnels for use in cytocentrifuges. Recommended for Thermo Scientific™ Shandon™ EZ Cytofunnel™ or other brands cytofunnels which have the same design.

	MARKING AREA	DESCRIPTION	P/N	
			CSW	SUPER GRADE
A	 White	Single printed circle Dia.: 13mm	0370-0003	0370-0103
	 White	Single printed circle Dia.: 13mm, coated	0370-0001	0370-0101
B	 White	Double printed circles Dia.: 13mm	0370-0004	0370-0104
	 White	Double printed circles Dia.: 13mm, coated	0370-0002	0370-0102
C	 White	Double printed circles Dia.: 13mm	0370-0020	0370-0120
	 White	Double printed circles Dia.: 13mm, coated	0370-0021	0370-0121

Microscope Slides with Printed Bracket

Specifically designed for use with the ThinPrep® and other brand processors.

	MARKING AREA	DESCRIPTION	P/N	
			CSW	SUPER GRADE
D	 White	Printed Bracket Dia.: 36mm, coated	0370-0013-16	0370-0113-16