

PCI Color-Plus™ PCI Microscope Slides **NEW**

Recommended for automatic IHC staining with the Roche Ventana IHC automated Stainer (privileged use), manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil, mark pens as well as the 2-D barcode based slide printers

Adhesion Microscope Slides for Printing Printers

Recommended for marking with the well-known brand printers, including the Inkjet IPS (Leica), Slide-Mate® (Thermo Fisher) and Signature® (Primera) slide printers

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophilic surface promotes the adhesion of specimens
- Color frosted area of approx. 20 mm depth, displays at one surface end
- Color frosted end is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requests and experiments



-50/PK
Packaging: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	Ground Edges 90° corners		Ground Edges 45° corners		Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7161-02	0303-7161-02	0312-7161-02	0302-7161-02	0314-7161-02	0304-7161-02
Green	0313-7161-08	0303-7161-08	0312-7161-08	0302-7161-08	0314-7161-08	0304-7161-08
Orange	0313-7161-10	0303-7161-10	0312-7161-10	0302-7161-10	0314-7161-10	0304-7161-10
Pink	0313-7161-12	0303-7161-12	0312-7161-12	0302-7161-12	0314-7161-12	0304-7161-12
White	0313-7161-16	0303-7161-16	0312-7161-16	0302-7161-16	0314-7161-16	0304-7161-16
Yellow	0313-7161-17	0303-7161-17	0312-7161-17	0302-7161-17	0314-7161-17	0304-7161-17

-72/PK
Packaging: 72pcs/box, 1,440pcs/carton

MARKING AREA	Ground Edges 90° corners		Ground Edges 45° corners		Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7171-02	0303-7171-02	0312-7171-02	0302-7171-02	0314-7171-02	0304-7171-02
Green	0313-7171-08	0303-7171-08	0312-7171-08	0302-7171-08	0314-7171-08	0304-7171-08
Orange	0313-7171-10	0303-7171-10	0312-7171-10	0302-7171-10	0314-7171-10	0304-7171-10
Pink	0313-7171-12	0303-7171-12	0312-7171-12	0302-7171-12	0314-7171-12	0304-7171-12
White	0313-7171-16	0303-7171-16	0312-7171-16	0302-7171-16	0314-7171-16	0304-7171-16
Yellow	0313-7171-17	0303-7171-17	0312-7171-17	0302-7171-17	0314-7171-17	0304-7171-17

PCI Color-Plus™ PCI Microscope Slides **NEW**


Recommended for automatic IHC staining with the Roche Ventana IHC automated Stainer (privileged use), manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections

*Markable with the conventional labeling system, pencil, mark pens as well as the 2-D barcode based slide printers

Adhesion Microscope Slides for Etching Printers


Recommended for marking with the slide etching labelers

- Dimensions: 25.0mm×75.0mm, Thickness: 1.1±0.1mm
- Pre-cleaned and hydrophilic surface promotes the adhesion of specimens
- Color frosted areas of approx. 20 mm depth, distribute on both sides but at the same slide end for marking
- Color frosted areas are resistant to all common chemicals and routine stains that are utilized in the laboratory
- A black color opposite the color side makes the marking clearer
- Beveled edges provide utmost safety during manual operations
- Various colors are supplied depending on your requests and experiments




-72/PK
Packaging: 72pcs/plastic box, 1,440pcs/carton


MARKING AREA	Ground Edges 90° corners		Ground Edges 45° corners		Beveled Edges 45° corners	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3271-02	0303-3271-02	0312-3271-02	0302-3271-02	0314-3271-02	0304-3271-02
Green	0313-3271-08	0303-3271-08	0312-3271-08	0302-3271-08	0314-3271-08	0304-3271-08
Orange	0313-3271-10	0303-3271-10	0312-3271-10	0302-3271-10	0314-3271-10	0304-3271-10
Pink	0313-3271-12	0303-3271-12	0312-3271-12	0302-3271-12	0314-3271-12	0304-3271-12
White	0313-3271-16	0303-3271-16	0312-3271-16	0302-3271-16	0314-3271-16	0304-3271-16
Yellow	0313-3271-17	0303-3271-17	0312-3271-17	0302-3271-17	0314-3271-17	0304-3271-17




Paraffin Embedded Sections




Frozen Sections




IHC



Human Blood Cells



Cervical Samples



Yeast Cells